



www.Urodry.com
1-888-391-7132

Warning: This product is made from latex which can cause allergic reactions in some individuals.

AlphaDry Instructions

ALPHADRY INCLUDES:

- ❖ AlphaDry Condom/Reservoir
- ❖ Velcro Strap
- ❖ Instructions/Care/Troubleshooting

Note: Use long-legged, brief-style underwear.

SEE ONLINE VIDEO
www.Urodry.com
Free help line
1-888-391-7132

Step 1



INITIAL PREPARATION:

Twist to open blue valve. Make sure penis and inside of condom are dry.

Step 2



ROLLING CONDOM DOWN:

Place thumbs over the inside bead and forefingers on the outside bead and begin to roll down like a condom.

Step 3



ROLL CONDOM DOWN TO INTERNAL VALVES:

Roll condom down until the internal valves are clearly showing.

Step 4



POSITION TO ROLL ALPHADRY ON PENIS:

Position condom at the head of penis as shown, with fingers of both hands on the bead and thumbs on the head of penis.

For Uncircumcised Men: Pull foreskin forward towards glans with one hand and position AlphaDry with other hand.

Option

Gel Strip can be used with the AlphaDry for added security, but always try the basic AlphaDry first.



Step 5



BEGIN TO ROLL CONDOM ONTO PENIS:

With your fingers, begin to roll the condom over the head of the penis.

Note: Make sure AlphaDry is the correct size. It should be snug when rolling it onto the penis. If it is too loose or too tight, you may order one AlphaDry replacement of a different size, at no extra charge, by calling 1-888-391-7132.

Step 6

UNROLL CONDOM TO BASE:

Unroll condom with both hands to base of penis.



More Instructions ➔

Step 7



WRAP ELASTIC STRAP BEHIND HEAD OF PENIS:

Hold the soft edge end of the strap about 1/2 inch away from the head of the penis angled slightly as shown. Begin to wrap the remainder of the strap right against the head of the penis with tension.

Step 8



LOCK STRAP DOWN:

After the strap is stretched, wrap strap on around and lock it down.

Step 9

IN PLACE ELASTIC STRAP

This is the correct position of the elastic strap.



Note: Positioning the strap just behind the head of the penis with tension is the key to preventing slipping and leaking. Adjust tension for comfort and to maintain dryness.

Option

COBAN wrap can be used in place of the Velcro strap, but always try the AlphaDry basic system first.



MAKE A HANDLE:

Pull a small amount of COBAN away from the roll. Then roll it on to itself so as to make a small roll. This will make it easy to remove later.



COBAN STRAP SHOWN IN PLACE:

This is the correct position of the COBAN strap.

Step 10



CLOSE BLUE VALVE

Close the blue valve on the AlphaDry and position straight down the leg in snug long-legged, boxer-type briefs.

Note: This position allows the one way valves to be free and open.

To purchase long-legged, boxer-briefs, visit www.Urodry.com or call us at 1-888-391-7132 for more information.

SEE ONLINE VIDEO
www.Urodry.com

Free help line
1-888-391-7132



You are now ready to resume life's normal activities with greater confidence.

— SEE WWW.URODRY.COM FOR: —

- ❖ Emptying the AlphaDry
- ❖ Troubleshooting
- ❖ Removing the AlphaDry
- ❖ Care and Helpful Hints