AlphaDry *Options*

The AlphaDry system provides you with options that help keep you dry and free again.

Compare and determine the options that are best for you. If you have questions or need more information to order, call toll free:

1-888-391-7132

Distributed by: www.UroDry.com

Warning: This product is made from latex which can cause allergic reactions in some individuals.




Standard

Note: Skin Shield[™] is available for skin protection and increased gripping.

STANDARD GRIPPING SYSTEM

AlphaDry use the natural gripping action of latex for the condom.

Option



OPTIONAL GEL STRIP

Applying the sticky Gel Strip behind head of the penis provides additional gripping power — and does not leave residue.



STANDARD SECURING


Velcro® elastic strap secures the AlphaDry, and can be reused several times.



OPTIONAL COBAN SECURING

The elastic 3M CobanTM self adherent wrap functions like tape but only sticks to itself.


Standard



STANDARD PLACEMENT

AlphaDry is placed somewhat horizontally inside the briefs.





Option



OPTIONAL PLACEMENT

The AlphaDry hangs straight down in snug, longer leg, boxer-style briefs.



OPTIONAL FLIP VALVE

Use the Flip Valve for the larger 2,000 ml reservoir, or for those who have trouble with manual dexterity.

